



Optimization
Makes Differences

OPM-293 CD05

Chemically Defined Cell Culture Medium

— For Biomanufacturing with HEK293 cell lines



OPM-293 CD05 is chemically defined, free of any animal-origin components, and contains no hydrolysates, proteins, growth factors or components of unknown composition. It can be used to reach & maintain high density suspension culture of HEK293 cell lines and to achieve highly efficient transient transfection. This liquid medium contains L-Glutamine.

Application

OPM-293 CD05 is intended for large scale manufacturing of therapeutic biomolecules, as well as for research purposes, but not for human or any therapeutic use.

Storage & Transportation

Store at 2~8°C, dark and dry
Ship at Room temperature (Liquid), Blue ice (Dry powder)

Shelf Life

OPM-293 CD05 Medium Liquid: 12 months

Quality Specifications

Specifications	OPM-293 CD05 Medium
Appearance	Red clear liquid
pH	7.0~7.4
Osmolality (mOsm/kg)	270~300
Solubility	—
Endotoxin (EUcells/mL)	<1.0
Sterility test	Negative

Cell Culture Conditions

37°C, 80% humidity, 5~8%CO₂
Shaker speed 110-150 rpm (amplitude: 50mm).

Cell Recovery

1. Rapidly thaw (<2 minute) a vial of frozen cells in a 37 °C water bath.
2. Transfer the entire contents aseptically into a 125mL shake flask containing 30 mL prewarmed OPM-293 CD05 Basal Medium.
3. Incubate at 37 °C in a humidified atmosphere of 5~8% CO₂ in air on a shaker (rotating at 115~125 rpm, amplitude: 50mm).

Cell Culture Passaging

1. Cultures should be passaged during the mid-log phase, approximately every 2~3 days. Be aware that different HEK293 cell lines may have different logarithmic growth phases which require individual calculation.

- Determine the viable cell density (VCD) and the cell viability.
- Proceed if $VCD = 3 \times 10^6 \sim 4 \times 10^6$ cells/mL & viability $\geq 95\%$. Otherwise troubleshoot conditions before continuing.
- Determine the correct volume of cell culture to inoculate a new flask at a starting cell density of $0.3 \times 10^6 \sim 0.6 \times 10^6$ cells/mL in prewarmed OPM-293 CD05 Basal Medium.

Medium Adaptation

Direct Medium Adaptation

- In most cases, HEK293 cell lines may be adapted directly from previous medium into OPM-293 CD05 Basal Medium.
- Adaptation should begin when cells are in mid-log phase, and viability $\geq 95\%$.
- Adaptation is completed when the cultures attain stable VCD of 3×10^6 cells/mL and viability $\geq 95\%$ within 3~4 days over at least 2~3 passages.

Sequential Medium Adaptation

- For certain HEK293 cell lines failing the direct medium adaptation, sequential adaptation is recommended.
- Passage the cells in current media for 2~3 passages to reach stable cell growth before beginning with medium adaptation.
- The adaptation instruction below provided below relies on maintaining the cell culture in mid-log growth phase by passaging the cells every 3 to 4 days. At least two passages at each adaptation step are recommended.
- Adaptation is completed when the cultures attain stable VCD of 3×10^6 cells/mL and viability $\geq 95\%$ within 3~4 days over at least 2~3 passages.

OPM-293 CD05 : current media (%)	Seeding density ($\times 10^6$ cells/mL)	Evaluation of cell growth	Acceptance criteria for next step
0:100	as usual	VCD & Viability	VCD $\geq 3 \times 10^6$ /mL, Viability $\geq 95\%$ over 2 passages
30:70	0.6	VCD & Viability	VCD $\geq 3 \times 10^6$ /mL, Viability $\geq 95\%$ over 2 passages
70:30	0.5	VCD & Viability	VCD $\geq 3 \times 10^6$ /mL, Viability $\geq 95\%$ over 2 passages
100:0	0.4	VCD & Viability	VCD $\geq 3 \times 10^6$ /mL, Viability $\geq 95\%$ over 2 passages

Recommended Transient Transfection Protocol

Recommended transfection condition

The optimal transfection conditions should be optimized case by case and may need to be established by DOE method. The following transfection conditions are only for reference:

VCD	$2 \times 10^6 \sim 4 \times 10^6$ cells/mL
DNA	0.7 ~ 2 mg/L
DNA/PEI ratio	1/2 ~ 1/6

Recommended transfection and feeding strategy

Feed	Feeding strategy
OPM-293 ProFeed	24h/48h post-transfection add 5% OPM-293 ProFeed

Recommended transfection procedure

Time line	Steps	Instruction
Preparation	1	Culture HEK293 cells used for transfection until stable VCD of $3 \times 10^6 \sim 4 \times 10^6$ cells/mL and viability $\geq 95\%$
Day -1	2	24 hours prior to transfection, split cells to a density of $3 \times 10^6 \sim 4 \times 10^6$ cells/mL and culture overnight in shake flasks.
Day 0	3	At the time of transfection, cell density should be $3 \times 10^6 \sim 4 \times 10^6$ cells/mL & with viability $\geq 95\%$. In case of higher VCD, dilute the cell suspension with fresh medium.
	4	Prepare the dilution of the expression vector plasmid DNA using OPM-293 CD05 Basal Medium. Mix carefully.
	5	Prepare a dilution of the PEI stock solution using OPM-293 CD05 Basal Medium. Mix carefully.
	6	Add diluted PEI to the diluted DNA plasmid solution, mix carefully. Incubate at room temperature 20 minutes.
	7	Add the DNA-PEI mixture to the cell suspension prepared in step 3 while swirling the culture to mix.
	8	Return the shake flask to 37°C incubator.
Day 1	9	16~24 hours post transfection, add 5% (v/v) OPM-293 ProFeed to the culture flask while slightly swirling. Return the flask to the 37°C incubator.
Day2-7	10	Maintain the glucose concentration above 4g/L. Harvest cells when the viability drops below 60%.



Optimization
Makes Differences

Order Information

Cell Culture Media

Name	Cat No.	Type	Volume	Description
OPM-293 CD05 Medium	81075-001	Liquid	1000mL	With L-Glutamine

High Performance Feeds

Name	Cat No.	Type	Volume	Description
OPM-293 ProFeed	F081918	Liquid	100mL	Protein-free feed
	F081918-001	Liquid	1000mL	

Other OPM-CD 293 Culture Medium

Name	Cat No.	Type	Volume	Description
OPM-293 CD03 DPM	91070-010	Dry powder	10L	Without L-Glutamine
	91070-050	Dry powder	50L	
OPM-293 CD03 Medium	81070-001	Liquid	1000mL	Without L-Glutamine
OPM-CD Trans293	P82019	Liquid	1000mL	With L-Glutamine

 Shanghai OPM Biosciences Co., Ltd.

OPM Headquarter: Building #28, 908 Ziping Road, Pudong, Shanghai, China
CDMO Development Center: Building #3, 100 Banxia Road, Pudong, Shanghai, China
Media & CDMO Manufacturing Center: C3&D3, No. 356, Zhengbo Road, Fengxian, China

+86 021-6818 2622
service@opmbiosciences.com
www.opmbio.com