

Sustained Viability Unlocks Extended Productivity Across CHO Transient Expression Hosts

Pargol Hashemi and Yunfen (Fiona) He
OPM Biosciences, Inc., Pleasanton, CA, USA

OPM
BIOSCIENCES

1. Introduction

Transient gene expression (TGE) in Chinese hamster ovary (CHO) cells is widely used for rapid production of recombinant proteins during early-stage biopharmaceutical research and development [1,2]. However, the performance of transient expression systems can depend strongly on the CHO host background, and commercially available platforms are often optimized for specific cell lines, such as CHO-S or CHO-K1-derived cells [1,2,3].

The OPM CD Trans® CHO Plus Expression System was developed as a streamlined, high-performance platform for transient protein production in CHO cells. Its FlashFection™ workflow enables direct addition of plasmid DNA and transfection reagent to the cell culture without separate dilution or pre-mixing steps, thereby reducing workflow complexity and simplifying transfection setup.

In this study, we evaluated the performance and host-cell compatibility of the OPM system in both CHO-S and CHO-K1-derived cells using two recombinant protein formats: an Fc-fusion protein and a monoclonal antibody (mAb). Performance was benchmarked against two leading commercially available transient expression systems optimized for the respective CHO host backgrounds.

Objective: Evaluate the productivity, culture longevity, and cross-host compatibility of the OPM CD Trans® CHO Plus Expression System in CHO-S and CHO-K1-derived cells relative to leading commercial transient expression systems.

2. Materials & Methods

Category	Description
Cell Lines	CHO-S and CHO-K1 derived cell lines
Expression Systems	OPM CD Trans CHO Plus Expression System, Competitor CHO-S-/CHO-K1-Optimized Expression Systems
Protein Models	Fc-based fusion protein (n=1: CHO-S; n=2: CHO-K1) and mAb (n=2: CHO-S & CHO-K1)
Seeding Density	CHO-S: 6×10^6 cells/mL; CHO-K1 derived : 9×10^6 (Competitor) and 6×10^6 cells/mL (OPM)
Culture Conditions	37 °C, 8% CO ₂ , orbital shaking; cells passaged 1 day prior to transfection
Transfection Conditions	Cell viability confirmed prior to use; plasmid DNA added according to manufacturer's protocol; DNA:Transfection Reagent ratio adjusted based on system
Analytical readouts	VCD (viable cell density), viability, titer, glucose (not shown here) and lactate

3. Simplified Protein Expression Workflow




Figure 1. Simplified workflow for the OPM CD Trans® CHO Plus Expression System. The FlashFection™ workflow enables direct addition of plasmid DNA and transfection reagent to the culture without separate dilution or pre-mixing steps, followed by scheduled Feed, Boost, Enhancer, and glucose supplementation.

4.1 Results: OPM vs. Competitor CHO-S Optimized System




Figure 2. OPM CD Trans® CHO Plus sustains viability and extends protein production in CHO-S cells. (A) Protein titers for Fc-fusion protein and mAb. (B) VCD and viability profiles throughout the culture. OPM maintained substantially higher viability through Day 10, supporting continued protein production and higher late-stage titers compared with the competitor CHO-S-optimized expression system.




Figure 2. Lactate profiles in CHO-S cells. OPM CD Trans® CHO Plus maintained low lactate levels throughout the culture for both the Fc-fusion protein and human IgG, while greater lactate accumulation was observed with the competitor CHO-S-optimized expression system.

4.2 Results: OPM vs. Competitor CHO-K1 Optimized System




Figure 3. OPM CD Trans® CHO Plus improves sustained productivity in CHO-K1-derived cells. (A) Protein titers and Day 9 normalized total protein for Fc-fusion protein and mAb. (B) VCD and viability profiles throughout the culture. (C) Lactate profiles. OPM maintained higher culture viability and sustained protein production through Day 9, resulting in substantially greater normalized total protein than the CHO-K1-optimized expression system. OPM also maintained low lactate levels, particularly during the later stages of culture, while greater lactate accumulation was observed with the competitor system.

5. Conclusions

- OPM CD Trans® CHO Plus demonstrated robust performance across both CHO-S and CHO-K1-derived host backgrounds and across two recombinant protein formats.
- Sustained culture viability extended the productive phase, enabling continued protein production and higher late-stage titers compared with host-optimized commercial expression systems.
- OPM demonstrated a substantial productivity advantage across both CHO host backgrounds and protein formats:

CHO Host	Protein	Comparison Day	OPM Fold Increase vs. Competitor (~)
CHO-S	Fusion protein	Harvest day	1.5 ×
	mAb		3.5 ×
CHO-K1-derived	Fusion protein		7 ×
	mAb		7 ×

- Together with the simplified FlashFection™ workflow, these results support OPM CD Trans® CHO Plus as a versatile platform for transient CHO protein expression.

Acknowledgements

We thank R&D and broader teams at OPM Biosciences (Shanghai & US) for their contributions to this experimental work and support.

References

- [1] Fu, Yushun et al. *Applied microbiology and biotechnology* vol. 108, 1 480. 4 Oct. 2024.
- [2] Reinhart, David et al. *Biotechnology journal* vol. 14, 3 (2019): e1700686.
- [3] Geisse, Sabine, and Cornelia Fux. *Methods in enzymology* vol. 463 (2009): 223-38.