



Maximize Your Protein Production with OPM-CHO Platform

OPM offers animal-origin-free basal media and feeds optimized for high productivity across different CHO cell lines including CHO-K1, CHOZN, CHO-DG44, and CHO-S. Our formulations support high cell densities, consistent titer improvements, and seamless GMP integration.

Why OPM-CHO Platform?

- Achieve an average 96% increase in titer while maintaining high cell viability
- Animal-origin-free and chemically defined for robust, reproducible results
- GMP-ready for seamless R&D-to-GMP transition
- Flexible packaging options (liquid bottles, buckets, or bags) and formats (DPM* or liquid)
- Robust intra- and inter-lot consistency (RSD <5%)

| Type | OPM-CHO Products | Available Sizes |
|------------------|--------------------|--|
| Basal Media | AltairCHO® | Liquid: 1000 mL DPM: 10/50/100 L |
| | SagiCHO™ | Liquid: 1000 mL DPM: 10/50/100 L |
| | StarCHO™ | Liquid: 1000 mL DPM: 10/50/100 L |
| | StarCHO™ Plus | Liquid: 500/1000 mL DPM: 10/50/500 L |
| | VegaCHO® | Liquid: 1000 mL DPM: 10/50/100 L |
| | HelixCHO™ | DPM: 10/50/100 L |
| Feeds | DenebCHO™ Feed | Liquid: 500/1000 mL DPM: 1/10/50 L |
| | SagiCHO™ Feed | Liquid: 1000 mL DPM: 10/50/100 L |
| | StarCHO™ Feed Plus | Liquid: 500/1000 mL DPM: 5/10/50 L |
| | HelixCHO™ Feed | Please inquire |
| | VectorCHO™ Feed | Please inquire |
| Feed Supplements | CDFS36™ Feed | Liquid: 500/1000 mL DPM: 1/10/50 L |
| | CDFS12™ Feed | Liquid: 500/1000 mL DPM: 5/10/50 L |

*DPM: Dry Powder Media

AltairCHO Basal Media and Feeds

Recommended for CHO-K1, CHOK1SV GS-KO, CHO-S, and CHO-DXB11




Figure 1. AltairCHO outperformed two global competitors in CHO-K1 mAb expression, showing a 31% titer increase with comparable VCD and viability. **(Left)** Titer comparison with two competitor products. **(Right)** VCD and viability profiles are shown.

SagiCHO Basal Media and Feeds

Recommended for CHO-K1, CHO-S, and CHO-GS




Figure 2. SagiCHO, combined with AltairCHO Feed, outperformed a commercial competitor in CHO-K1 mAb expression — showing a 130% titer increase with higher late-stage VCD and viability. **(Left)** Titer comparison with competitor product. **(Right)** VCD and viability profiles are shown.

StarCHO Basal Media and Feeds

Recommended for CHOZN[®] GS^{-/-}




Figure 3. StarCHO outperformed a commercial competitor in CHOZN GS^{-/-} mAb expression, showing a 130% titer increase with comparable VCD and viability. **(Left)** Titer comparison with competitor product. **(Right)** VCD and viability profiles are shown.

VegaCHO Basal Media and Feeds

Recommended for CHO-DG44





Figure 4. VegaCHO outperformed 2 commercial competitors in mAb expression, showing a 81% titer increase with comparable VCD and viability. **(Left)** Comparison of titer with two competitor products. **(Right)** VCD and viability profiles are shown.

High Titers Powered by OPM-CHO Products

OPM has completed 35 protein titer improvement projects using OPM-CHO products, achieving an average **96%** increase — and up to **650%** in a single project.



Media Selection Guide

Select the optimal media and feed for your cell line.

| Cell Lines | Basal Medium | Feeds | Feed Supplement |
|---------------------------------|-------------------------|--------------------------------|-----------------|
| CHO-K1 (ATCC) | AltairCHO SagiCHO | DenebCHO Feed SagiCHO Feed | CDFS36 |
| CHO-K1 (ECACC) | AltairCHO | DenebCHO Feed SagiCHO Feed | CDFS36 |
| CHOK1SV GS-KO (Lonza) | AltairCHO | DenebCHO Feed | CDFS36 |
| CHOZN GS ^{-/-} (Merck) | StarCHO StarCHO Plus | StarCHO Feed Plus | CDFS12 |
| CHOZN CHO-K1 (Merck) | HelixCHO | HelixCHO Feed | CDFS36 |
| Horizon CHO-GS | SagiCHO | SagiCHO Feed | CDFS36 |
| CHO-DG44 | VegaCHO | VectorCHO Feed | CDFS36 |
| CHO-S | SagiCHO AltairCHO | SagiCHO Feed | CDFS36 |
| CHO-DXB11 | AltairCHO | SagiCHO Feed VectorCHO Feed | CDFS36 |